



# GP Detergent

## All-Purpose Manual Dish Wash Detergent

5 Litres	- P563
20 Litres	- P566
200 Litres	- P569

### DESCRIPTION:

GP Detergent is a general purpose detergent, suitable for cleaning all types of kitchen equipment and surfaces. It is excellent on plates, cutlery, crockery and greasy pots & pans.

It is phosphate-free and ideal for cleaning food and meat preparation appliances and utensils. It is biodegradable.

### DIRECTIONS:

Dilute GP Detergent 1 part to 150 parts water. Add additional detergent for greasy dishes. For greater economy, use a tap proportioner to dispense product. Rinse with potable water after use.

**CERTIFICATION:** NZFSA Approved C31 (All animal product except dairy)

### PHYSICAL AND CHEMICAL PROPERTIES:

<b>Appearance:</b>	Clear Green
<b>Odour:</b>	Odourless
<b>Specific Gravity:</b>	1.00
<b>pH:</b>	7.5 +/- 0.5

### FIRST AID MEASURES:

#### Ingestion:

Immediately rinse mouth with water. If swallowed do not induce vomiting. Give water to drink. Never give anything by mouth to an unconscious person. Seek immediate medical attention.

#### Eye Contact:

Immediately flush eyes with large amounts of water for at least 15 minutes while holding eyelids open. Transport to the nearest medical facility for additional treatment. Do not rub eyes or keep eyes closed.

#### Skin Contact:

Remove contaminated clothing. Flush exposed area with water and follow by washing with soap if available.

#### Inhalation:

Remove the affected person out to a ventilated area. If not breathing, give artificial respiration. If breathing is difficult, give oxygen. Seek medical attention.